

REALQUALITY RQ-HHV 8 STANDARD	
CE IVD	Ready-to-use quantified standards for quantification of <i>Human Herpes virus type 8 (HHV 8)</i>
Code RQ-18	
Description	The REALQUALITY RQ-HHV 8 STANDARD is to be used exclusively in association with the REALQUALITY RQ-HHV 8 kit, code RQ-17, for quantitative testing. The standard consists of DNA fragments containing the gene region of interest and therefore is not dangerous to the user.
CND/EDMS code	W0105080904/14500190
Reagent status	Ready-to-use
Packaging	4 x 60 µL 4 calibration points are provided, each sub-aliquoted in 3 x 20 µL, in sufficient quantity to do 12 standard curves.
Stability	12 months Each standard aliquote can undergo a maximum number of 3 freeze/thaw cycles: further cycles compromise the stability of the product and may alter its concentration.
Storage	-30°C/-20°C It is recommended to freeze the standard after each use.
Concentration	STANDARD 1: 20 viral genome copies /µL STANDARD 2: 200 viral genome copies /µL STANDARD 3: 2.000 viral genome copies /µL STANDARD 4: 20.000 viral genome copies /µL
Instructions for use	To obtain a standard curve in the range of 10 ² -10 ⁵ copies of DNA, use 5 µL for reaction, strictly adhering to the indications given in the user manual of the REALQUALITY RQ-HHV 8 kit, code RQ-17. It is recommended to vortex and centrifuge each tube before use.

ORDER INFORMATION		
Code	Product	PKG
RQ-18-SM	REALQUALITY RQ-HHV 8 STANDARD	4 x 60 µL
In combination with the product:		
RQ-17-4M	REALQUALITY RQ-HHV 8	50 tests
RQ-17-6M		100 tests